


Designed exclusively for animals  
in pursuit of easy to use



# Tierrett Long Foley Catheter

Registration No. of Medical Equipment for Animal Use: 26 Douyaku 3980



**Recommend to large male dogs**  
For dogs whose urine leaks at 8Fr!

## Lineup

Size	O.D.	Balloon capacity	Effective length	Product cord	JAN cord
6Fr	2.0mm	1mL	650mm	800-001-3467	4542187191356
			850mm	800-001-3469	4542187191370
8Fr	2.7mm	3mL	650mm	800-001-3470	4542187191387
			850mm	800-001-3472	4542187191400
10Fr	3.3mm	3mL	650mm	800-001-6679	4542187223446
			850mm	800-001-6680	4542187223453

**New!!**

5-pack box/EOG sterilized

⚠ Please note that the product specifications are subject to change without notice to reflect the results of our consistent R&D efforts.  
⚠ Read the package insert before using the product.



**CREATE MEDIC CO., LTD.**





Made of human-body-friendly silicon



Moderately elastic hydrophilic stylet for improve insertability



Easy removal

Easy-to-handle lightweight /compact joint



Reduced weight to 50% of conventional products and 55% of length



Multi-cap that can be used in 3 ways



Catheter tip

Lure tip

Stopcock



**CREATE MEDIC CO., LTD.**

<https://www.createmedic.co.jp/english/>

【Contact information】

Animal Medical Business Department E-mail : oem-infoch@createmedic.co.jp

【Head office】 2-5-25 Chigasaki Minami, Tsuzuki-ku, Yokohama, Kanagawa, Japan



Control No : 30-1826-20201014