



1. The design of catheter case.

Since the catheter case has not only good performance as medical device but also fashionable design, it does not look like medical device at first glance.



3. The strap with hook.

Since this strap has hook, it can be hooked anywhere and easy to get it out.

5. Reusable type.

After using the catheter, if it is washed and stored in the sterile liquid it can be used for one (1) month as reusable product.

2. High flow catheter.

This catheter made from all silicone that has a good performance of biocompability. Since its lumen is wider than the other Self-Catheterization products, the time of urethral catheterization can be shortened.



4. The case cap.

Since the case cap is designed as an oval to be held easily, woman can open and fasten it easily. When the user intends to open, since twist the cap and it will be come up, the catheter can be picked up easily. On the other hand, when it is closed, only push the cap and it comes to fix the cap and its case. Moreover, the marks where are on the cap and case indicate correct position to be fastened firmly.



Specification

Product number	Size	Outer diameter	Length	Side hole	Contents
800 000 9913	9Fr	3.0mm	125mm	2 holes	Catheter × 1piece Catheter case × 1piece Strap × 1piece
800 000 9915	12Fr	4.0mm	125mm		
800 000 9916	14Fr	4.7mm	125mm		
800 001 1323	12Fr (Hard Type)	4.0mm	125mm		

Package: 5 sets per box, Non-sterile.

Dr. Takayuki Sugimura (Urologist) and Ms. Junko Tanaka (Nurse and Urinate Function Coordinator), Urology Department of St Luke's International Hospital, are acknowledged for their generous offer of direction and advices in development of the device.

Please note that a partial amendment in specifications is made occasionally, as studies and developments been performed continuously.

Be sure to read through the instruction for use attached to the device before



005